


Piedro reverse lasted correction footwear has the purpose to correct feet with a too large varus and/or adduction position. This correction takes place due to the special last shape. Inside, the shoe is equipped with a special insole that in combination with a medial extended stiffener which forms a curved wall through the connection with the toe cap avoids that the foot comes in a further adduction position. At the same time the abduction of the feet will be stimulated. This product is especially very important for children who are born with clubfeet; treatment with the reverse lasted shoes will help them to develop their feet.

Reverse lasted footwear

The Piedro reverse lasted footwear is recommended to be used by:

- > Corrects the peak position in the upper tarsal joint (Pes equino-varus);
- > Corrects the varus tilting in the bottom tarsal joint and the further metatarsal area (Pes varus);
- > Stimulates the abduction of the mid- and forefoot;
- > Stabalises the post-operative clubfeet;
- > Improves the toe-off of the foot.




This is our Piedro Reverse lasted stock brochure. For additional possibilities see Piedro Select.

Piedro Verkooporganisatie B.V.
Jan van der Heijdenstraat 18-20
NL-2665 JA Bleiswijk

T +31 10 521 52 11
F +31 10 521 49 61
E sales@piedro.com
www.piedro.com

> Piedro. Always one step ahead



>Stock collection



>Flexibility

>Stability

>Durable

>Comfortable

Productrange
Piedro Reverse lasted footwear


>Stock collection



→18 - 31 5 

Art no. 2912.5600




→18 - 33 5 

Art no. 2936.5600

All models are available directly from stock. Do you wish to make your own shoe? Please look in our Piedro Select brochure.

>Reverse lasted footwear



→18 - 33 5 

Art no. 2978.6509



→18 - 33 5 

Art no. 2987.5600



→18 - 33 5 

Art no. 2975.1602

>Sizing Chart

Piedro, Cemented making Reverse lasted (heel/ball correction inlay)			
Size	Inner length in mm.	Angle	Width 5 in cm.
20	133	11,50°	14,85
21	140	11,50°	15,25
22	146	11,00°	15,65
23	153	11,00°	16,05
24	160	10,50°	16,45
25	166	10,50°	16,85
26	173	10,00°	17,25
27	179	10,00°	17,65
28	186	10,00°	18,05
29	193	9,50°	18,45
30	199	9,50°	18,80
31	206	9,00°	19,15
32	213	9,00°	19,50
33	219	8,50°	19,90
34	226	8,50°	20,30
35	233	8,50°	20,70

Piedro, Stitched down making Reverse lasted (heel/ball correction inlay)			
Size	Inner length in mm.	Angle	Width 5 in cm.
18	117	11,00°	14,30
19	123	11,00°	14,70
20	130	10,50°	15,10
21	137	10,50°	15,50
22	143	10,50°	15,90
23	150	10,00°	16,30
24	157	10,00°	16,70
25	163	10,00°	17,10
26	170	9,50°	17,50
27	177	9,50°	17,90
28	183	9,50°	18,30
29	190	9,00°	18,70
30	197	9,00°	19,10
31	203	9,00°	19,50
32	210	8,50°	19,90
33	217	8,50°	20,30

If you wish insole patterns, please call us (+31) (0) 10 521 5211, or send an e-mail to sales@piedro.com

